

# Standardization of Terminology of Lower Urinary Tract Function: Pressure-Flow Studies of Voiding, Urethral Resistance, and Urethral Obstruction

Derek Griffiths, Klaus Höfner, Ron van Mastrigt, Harm Jan Rollema, Anders Spångberg, and Donald Gleason

*International Continence Society Subcommittee on Standardization of Terminology of Pressure-Flow Studies, coordinating chairman Anders Mattiasson*

## LIST OF CONTENTS

1. Introduction
  2. Evaluation of micturition
    - 2.1 Pressure-flow studies
      - 2.1.1 Pressure and flow rate parameters
      - 2.1.2 Flow delay
      - 2.1.3 Presentation of results
    - 2.2 Urethral resistance and bladder outlet obstruction
      - 2.2.1 Urethral function of during voiding
      - 2.2.2 Urethral resistance
      - 2.2.3 Urethral activity
      - 2.2.4 Bladder outlet obstruction
      - 2.2.5 Methods of analyzing pressure-flow plots
    - 2.3 Detrusor contractility during micturition
  3. Additional symbols
- Appendix: ICS standard for digital exchange of pressure-flow study data

**Key words:** terminology; pressure flow analysis; urethral resistance; urethral obstruction; URR; ICS standards

## 1. INTRODUCTION

This report has been produced at the request of the International Continence Society. It was approved at the twenty-fifth annual meeting of the Society in Sydney, Australia.

The 1988 version of the collated reports on standardization of terminology, which appeared in *Neurourology and Urodynamics* 7:403–427, contains material relevant to pressure-flow studies in many different sections. This report is a revision and expansion of sections 4.2 and 4.3 and parts of sections 6.2 and 7 of the 1988




Fig. 1. Diagram of a pressure-flow study with nomenclature recommended in this report.

report. It contains a recommendation for a provisional standard method for defining obstruction on the basis of pressure-flow data.

## 2. EVALUATION OF MICTURITION

### 2.1 Pressure-flow studies

At present, the best method of analyzing voiding function quantitatively is the pressure-flow study of micturition, with simultaneous recording of abdominal, intravesical, and detrusor pressures and flow rate (Fig. 1).

Direct inspection of the raw pressure and flow data before, during, and at the end of micturition is essential because it allows artifacts and untrustworthy data to be recognized and eliminated. More detailed analyses of pressure-flow relationships, described below, are advisable to aid diagnosis and to quantify data for research studies.

The flow pattern in a pressure-flow study should be representative of free-flow studies in the same patient. It is important to eliminate artifacts and unrepresentative studies before applying more detailed analyses.

Pressure-flow studies contain information about the behavior of the urethra and the behavior of the detrusor. Section 2.2 deals with the urethra. Detrusor function is considered in section 2.3.

### 2.1.1 Pressure and flow rate parameters

**Definitions.** See Fig. 1 and Table II; see also Table II in the 1988 version of the collated standardization reports.

*Maximum flow rate* is the maximum measured value of the flow rate. Symbol:  $Q_{\max}$ .

*Maximum pressure* is the maximum value of the pressure measured during a pressure-flow study. Note that this may be attained at a moment when the flow rate is zero. Symbols:  $p_{\text{abd. max}}$ ,  $p_{\text{ves. max}}$ ,  $p_{\text{det. max}}$ .

*Pressure at maximum flow* is the pressure recorded at maximum measured flow rate. If the same maximum value is attained more than once or if it is sustained for a period of time, then the point of maximum flow is taken to be where the detrusor pressure has its lowest value for this flow rate; abdominal, intravesical, and detrusor pressures at maximum flow are all read at the same point. Flow delay (see 2.1.2) may have a significant influence and should be considered. Symbols:  $p_{\text{abd. Qmax}}$ ,  $p_{\text{ves. Qmax}}$ ,  $p_{\text{det. Qmax}}$ .

*Opening pressure* is the pressure recorded at the onset of measured flow. Flow delay should be considered. Symbols:  $p_{\text{abd. open}}$ ,  $p_{\text{ves. open}}$ ,  $p_{\text{det. open}}$ .

*Closing pressure* is the pressure recorded at the end of measured flow. Flow delay should be considered. Symbols:  $p_{\text{abd. clos}}$ ,  $p_{\text{ves. clos}}$ ,  $p_{\text{det. clos}}$ .

*Minimum voiding pressure* is the minimum pressure during measurable flow (see Fig. 1). It may be, but is not necessarily, equal to the opening pressure or the closing pressure. For example, the symbol for minimum voiding detrusor pressure is  $p_{\text{det. min. void}}$ .

**2.1.2 Flow delay.** When a pressure-flow study is performed, the flow rate is measured at a location downstream of the bladder pressure measurement, so the flow rate measurement is delayed. The delay is partly physiological, but it also depends on the equipment. It may depend on the flow rate.

When considering pressure-flow relationships, it may be important to take this delay into account, especially if there are rapid changes in pressure and flow rate. In current practice an average value is estimated by each investigator from observations of the delay between corresponding pressure and flow rate changes in a number of actual studies. Values from 0.5 to 1.0 s are typical.

**Definition.** *Flow delay* is the time delay between a change in bladder pressure and the corresponding change in measured flow rate.

**2.1.3 Presentation of results.** Pressure-flow plots and the nomograms used for analysis should be presented with the flow rate plotted along the X-axis and the detrusor pressure along the Y-axis (see Fig. 2).

**Specify.** The value of the flow delay that is used.




Fig. 2. Plot of detrusor pressure against flow rate during voiding (broken curve), providing an indication of the urethral resistance relation (URR). The continuous smooth curve is an estimate of the passive urethral resistance relation.

## 2.2 Urethral resistance and bladder outlet obstruction

**2.2.1 Urethral function during voiding.** During voiding, urethral function may be

- a) Normal;
- b) Obstructive, as a result of:
  - i) overactivity;
  - ii) abnormal structure.

Obstruction due to urethral overactivity occurs when the urethral closure mechanism contracts involuntarily or fails to relax during attempted micturition in spite of an ongoing detrusor contraction. Obstruction due to abnormal structure has an anatomical basis, e.g., urethral stricture or prostatic enlargement.

**2.2.2 Urethral resistance.** Urethral resistance is represented by a relation between pressure and flow rate, describing the pressure required to propel any given flow rate through the urethra. The relation is called the *urethral resistance relation* (URR).

An indication of the urethral resistance relation is obtained by plotting detrusor pressure against flow rate. The most accurate procedure, which requires a computer or an x/y recorder, is a quasi-continuous plot showing many pairs of corresponding pressure and flow rate values (Fig. 2). A simpler procedure, which can be performed by hand, is to plot only 2 or 3 pressure-flow points connected by straight lines; for example, the points of minimum voiding pressure and of maximum flow may be selected. No matter how the plot is made, flow delay should be considered.

A further simplification is to plot just one point showing the maximum flow rate and the detrusor pressure at maximum flow. Flow delay should be considered.

Methods of analyzing pressure-flow plots are further discussed below.

**2.2.3 Urethral activity.** Ideally, the urethra is fully relaxed during voiding. The urethral resistance is then at its lowest and the detrusor pressure has its lowest value for any given flow rate. Under these circumstances the urethral resistance relation is defined by the inherent mechanical and morphological properties of the urethra and is called the *passive urethral resistance relation* (Fig. 2).

TABLE I. Methods of Analyzing Pressure-Flow Plots

Method	Aim	Number of $p/Q$ points	Assumed shape of URR	Number of parameters	Number of classes or continuous
Abrams/Griffiths nomogram [1]	diagnosis	1	n/a	n/a	3
Spångberg nomogram [10]	diagnosis	1	n/a	n/a	3
URA [2,6]	resistance	1	curved	1	continuous
linPURR [8]	resistance	1*	linear	1	7
Schäfer PURR [7]	resistance	many	curved	2	continuous
CHESS [3]	resistance	many	curved	2	16
OBI [4]	resistance	many	linear	1	continuous
Spångberg et al. [10]	resistance	many	linear or curved	3	continuous + 3 categories
DAMPF [9]	resistance	2	curved	1	continuous
A/G number [5]	resistance	1	linear	1	continuous

\*Schäfer uses 2 points to draw a linear relation, but the point at maximum flow determines the resistance grade.

Urethral activity can only increase the detrusor pressure above the value defined by the passive urethral resistance relation. Therefore, any deviations of the pressure-flow plot from the passive urethral resistance relation toward higher pressures are considered to be a result of activity of the urethral or periurethral muscles, striated or smooth.

**2.2.4 Bladder outlet obstruction.** Obstruction is a physical concept that is assessed from measurements of pressure and flow rate made during voiding. Whether due to urethral overactivity or to abnormal structure, obstruction implies that the urethral resistance to flow is abnormally elevated. Because of natural variation from subject to subject, there cannot be a sharp boundary between normal and abnormal. Therefore, the definition of abnormality requires further elaboration.

**2.2.5 Methods of analyzing pressure-flow plots.** The results of pressure-flow studies may be used for various purposes—for example, for objective diagnosis of urethral obstruction or for statistical testing of differences in urethral resistance between groups of patients. For these purposes, methods have been developed to quantify pressure-flow plots in terms of one or more numerical parameters. The parameters are based on aspects such as the position, slope, and curvature of the plot. Some of these methods are intended primarily for use in adult males with possible prostatic hypertrophy.

Some methods of analysis are shown in Table I.

*Quantification of urethral resistance.* In all current methods, urethral resistance is derived from the relationship between pressure and flow rate. A commonly used method of demonstrating this relationship is the pressure-flow plot. The lower pressure part of this plot is taken to represent the passive urethral resistance relation (see Fig. 2). In general, the higher the pressure for a given flow rate, and/or the steeper or more sharply curved upward this part of the plot, the higher the urethral resistance. The various methods differ in how position, slope, and/or curvature of the plot are quantified and how and whether they are combined. Some methods grade urethral resistance on a continuous scale; others grade it in a small number of classes (Table I). If there are few classes, small changes in resistance may not be detected. Conversely, a small change on a continuous scale may not be clinically relevant.

Some methods result in a single parameter; others result in 2 or more parameters (Table I). A single parameter makes it easy to compare different measurements. A larger number of parameters makes comparison more difficult but potentially gives higher accuracy and validity. If there are too many parameters, however, accuracy may be compromised by poor reproducibility.

*Choice of method.* Some methods in Table I are intended primarily to quantify urethral resistance. Others are intended only for the diagnosis of obstruction. Methods that quantify urethral resistance on a scale can also be used to aid diagnosis of obstruction by comparison with cut-off values. In every case an equivocal zone may be included.

Because of their underlying similarity, all the above methods classify clearly obstructed and clearly unobstructed pressure-flow studies consistently, but there is some lack of agreement in a minority of cases with intermediate urethral resistance.

One method of analyzing pressure-flow studies may be more useful for a particular purpose. In selecting a method, investigators should consider carefully what their aims are and which method is best suited to their purpose.

*Identification of optimum methods.* For a subsequent report, the International Continence Society will compare the above methods with each other and may also develop new methods, with the aim of reaching a consensus on their use. The Society will continue to seek better ways of clinically validating these methods. The following procedure has been chosen.

Making use of good-quality data stored in digital format, the following databases will be examined:

1. Pressure-flow studies in untreated men with lower urinary tract symptoms and signs suggestive of benign prostatic obstruction;
2. Pressure-flow studies repeated after a time interval with no intervention;
3. Pressure-flow studies before and after TURP;
4. Pressure-flow studies before and after alternative therapeutical intervention that causes a small change in urethral resistance.

Database 1 will be used to determine which existing or new methods adequately describe the actual pressure-flow plots of male patients with lower urinary tract symptoms. Database 2 will be used to determine the reproducibility of the various methods. Database 3 will be used to determine in which groups of patients TURP significantly reduces urethral resistance, and hence which patients are indeed obstructed. Database 4 will be used to test the various methods' sensitivity to small changes of urethral resistance.

On the basis of these analyses, the International Continence Society will attempt to identify

- a) A simple and reproducible method with high validity of diagnosing obstruction;
- b) A sensitive and reproducible method with high validity of measuring urethral resistance and changes in resistance.

*Provisional recommendation.* Pending the results of these procedures, it is recommended that investigators reporting pressure-flow studies in adult males, particularly subjects with benign prostatic hyperplasia, use one simple standard method of analysis in addition to any other method that they have selected, so that results




Fig. 3. Provisional ICS method for definition of obstruction. The points represent schematically the values of maximum flow rate and detrusor pressure at maximum flow for 9 different voids, 3 in each class.

from different centers can be compared. For this provisional method it is recommended that urethral resistance be specified by the maximum flow rate and the detrusor pressure at maximum flow, i.e., by the pair of values ( $Q_{\max}$ ,  $p_{\text{det.Qmax}}$ ). A provisional diagnostic classification may be derived from these values as follows:

- If  $(p_{\text{det.Qmax}} - 2 \cdot Q_{\max}) > 40$ , the pressure-flow study is obstructed;
- If  $(p_{\text{det.Qmax}} - 2 \cdot Q_{\max}) < 20$ , the pressure-flow study is unobstructed;
- Otherwise, the study is equivocal.

In these formulas, pressure and flow rate are expressed in cm H<sub>2</sub>O and ml/s respectively. This method is illustrated graphically in Figure 3. It may be referred to as the *provisional ICS method for definition of obstruction*.

The equivocal zone of the provisional method (Fig. 3) is similar but not identical to those of the Abrams-Griffiths and Spångberg nomograms and to the region defining linPURR grade II. For micturations with low to moderate flow rates it is consistent with cut-off values used to define obstruction in the URA and CHES methods.

### 2.3 Detrusor contractility during micturition

During micturition the detrusor may be:

- acontractile;
- underactive;
- normal.

An *acontractile detrusor* is one that cannot be demonstrated to contract during urodynamic studies.

An *underactive detrusor* produces a contraction of inadequate magnitude and/or duration to effect complete bladder emptying in the absence of urethral obstruction. (Concerning the elderly, see below.) Both magnitude and duration should be considered in the evaluation of detrusor contractility.

A *normal detrusor* in the absence of obstruction, produces a contraction that will effect complete bladder emptying. Detrusor contractility in the elderly may need special consideration.

For a given detrusor contraction, the magnitude of the recorded pressure rise

will depend on the outlet resistance. In general, the higher the detrusor pressure and/or the higher the flow rate, the stronger the detrusor contraction. The magnitude of the detrusor contraction may be quantified approximately by means of a nomogram applied to the pressure-flow plot or by calculation.

### 3. ADDITIONAL SYMBOLS

Qualifiers that can be used to form symbols for variables relevant to voiding are shown in Table II. These are additions to those in Table II of the 1988 standardization report.

**TABLE II. Qualifiers that Can Be Used to Indicate Pressure and Flow Variables Relevant to Voiding\***

Qualifiers	
At maximum flow	Qmax
During voiding	void
Opening	open
Closing	clos
Examples	
$P_{det.Qmax}$	= detrusor pressure at maximum flow
$P_{det.min.void}$	= minimum voiding detrusor pressure
$P_{ves.open}$	= intravesical opening pressure
$P_{ves.clos}$	= intravesical closing pressure

\*When possible, qualifiers should be printed as subscripts (see above). Note that the preferred symbol for pressure is lower-case  $p$ , while the symbol for flow rate is uppercase  $Q$ .

### REFERENCES

- Abrams PH, Griffiths DJ (1979): The assessment of prostatic obstruction from urodynamic measurements and from residual urine. *British Journal of Urology* 51:129–134.
- Griffiths D, Van Mastrigt R, Bosch R (1989): Quantification of urethral resistance and bladder function during voiding, with special reference to the effects of prostate size reduction on urethral obstruction due to benign prostatic hypertrophy. *Neurourol Urodynam* 8:17–27.
- Höfner K, Kramer AEJL, Tan HK, Krah H, Jonas U (1995): CHES classification of bladder outflow obstruction: a consequence in the discussion of current concepts. *World J Urol* 13:59–64.
- Kranse M, Van Mastrigt R (1991): The derivation of an obstruction index from a three-parameter model fitted to the lowest part of the pressure flow plot. *J Urol* 145:261A.
- Lim CS, Abrams P (1995): The Abrams-Griffiths nomogram. *World J Urol* 13:34–39.
- Rollema HJ, Van Mastrigt R (1992): Improved indication and follow-up in transurethral resection of the prostate (TUR) using the computer program CLIM. *J Urol* 148:111–116.
- Schäfer W (1983): The contribution of the bladder outlet to the relation between pressure and flow rate during micturition. In Hinman F Jr (ed), "Benign Prostatic Hypertrophy." New York: Springer Verlag, pp. 470–496.
- Schäfer W (1990): Basic principles and clinical application of advanced analysis of bladder voiding function. *Urol Clin N Am* 17:553–566.
- Schäfer W (1995): Analysis of bladder-outlet function with the linearized passive urethral resistance relation, linPURR, and a disease-specific approach for grading obstruction: from complex to simple. *World J Urol* 13:47–58.
- Spångberg A, Teriö H, Ask P, Engberg A (1991): Pressure/flow studies preoperatively and postoperatively in patients with benign prostatic hypertrophy: estimation of the urethral pressure/flow relation and urethral elasticity. *Neurourol Urodyn* 10:139–167.

**Appendix: ICS Standard for Digital Exchange of Pressure-Flow Study Data\*****LIST OF CONTENTS**

- A1. Introduction
- A2. General description of signal storage
- A3. Variable values and types
- A4. General structure of file and records
- A5. Definitions of record types
- A6. Signals and information to be stored: minimal specification and optional extensions
- A7. Typical file structure
- A8. Acknowledgments
- A9. Addendum: signal ID's

**A1. INTRODUCTION**

To facilitate exchange of digital urodynamic data a standard file format is required. In this document an ICS standard is summarized. Its primary purpose is to enable data from pressure-flow studies to be exchanged. Enough detail is given to allow exchange of other urodynamic data. Extensions may be made as described in section A5.

**A2. GENERAL DESCRIPTION OF SIGNAL STORAGE**

For each pressure-flow study, urodynamic signals sampled equidistantly with an A/D converter, and other associated information, are stored in one binary MS-DOS compatible file on a 5.25- or 3.5-in floppy disc. The stored signals start 10 s before a detectable change in the flow rate signal and continue until 10 s after the flow rate has finally returned to the baseline. Whenever possible, all signals should be stored at the same sample rate. In this case, all signals have the same length, i.e., they contain the same number of bytes.

The file name and extension are ABCDEFGH.ICS, where ABCDEFGH stands for a unique measurement identification string. In the case of a multicenter study requiring exchange of data, "unique" implies that this string is defined by the coordinating center. In the ICS-BPH study, for instance, this would be a number consisting of (first) 3 digits for center number, (next) 3 digits for patient number, and (finally) 2 digits for identification of successive measurements made on this one patient in one or more sessions.

The file consists of a number of records of various types. In the ICSMFF proposal on which this document is partly based (see Acknowledgments), 14 different types of records are defined, numbered -1 and 1 to 13. Existing record types may not be modified, but extensions may be implemented by defining new types of record. For ICS purposes, 6 additional types of record, types 14 to 19, are required, and these

are defined in section A5. In addition, the structures of records of types 8 and 9 are further elaborated in section A5.

### A3. VARIABLE VALUES AND TYPES

In this report, actual values for bytes or words may be given in hexadecimal notation as follows: Hex:DD for 8-bit bytes, Hex:DDDD for 16-bit words, and so on, where D stands for a hexadecimal digit.

The following variable types are used in the definitions in the following sections:

byte	:unsigned 8-bit value
word	:unsigned 16-bit value
integer	:signed 16-bit value
dword	:unsigned 32-bit value
string[N]	:N bytes including a terminating Hex:00
word [N]	:N words

### A4. GENERAL STRUCTURE OF FILE AND RECORDS

Each file ABCDEFGH.ICS consists of a number of records. The number of records is not predefined. Records may be of different types, containing different kinds of data. For example, a record may contain a description of a urodynamic signal or information about the patient.

Each record starts with a descriptor containing the record type, the size of the data block in bytes, and a checksum. The descriptor is followed by the block containing the actual data (see Figure 4).

The record type and the block size are integers. The checksum is an integer such that the 16-bit sum of the record type, the data block size, and the checksum is zero.

### A5. DEFINITIONS OF RECORD TYPES

The record type defines what kind of data is contained by the data block in the record. The only constraint is that the record type is unique.

Backward compatibility is preserved by the requirement that previously defined record types may not be modified. Extensions are implemented by defining new types of record. New types of signal may be implemented by defining new signal ID's, together with signal names and units. Forward compatibility will be ensured if a user does not try to handle unknown types of record or unknown signal ID's.

A register of the types of record in use, and a register of signal ID's, signal names, and signal units, will be administered centrally by the ICS. This standard contains initial versions of the registers. Proposed extensions should be communicated to the ICS for registration.

#### Default Values

Empty fields are not permitted. Occasionally, the value for a certain field may not be available. In this case a default value should be stored (see below).




Fig. 4. Schematic structure of file and records.

**Signal ID (See A9: Addendum)**

The following sections describe the types of record needed to implement this ICS standard. Other types are described in the document ICS Measure File Format (ICSMFF) version 1.00, dated 14.10.93 (see Acknowledgments).

**RECORD TYPE 8: THE MEASUREMENT MARKER RECORD**

This type of record defines a marker within the measurement and an associated comment.

Record type:	8	:integer
Blocksize:	1 + N + 1 + 4 + M	:integer
Checksum:	-(8 + blocksize)	:integer
Data block:	marker type	:byte
	name	:string[N]
	signal id	:byte
	position	:dword
	comment	:string[M]

## 12 Griffiths et al.

The following information should be stored in these fields:  
In the field:

marker type	0 = default value (other types are defined in the ICSMFF document)
name	empty string as default value (empty string contains only the terminating Hex:00, so that $N = 1$ )
signal id	as in A9: Addendum
position	the number of the sample (16-bit word) with which the marker and comment are associated
comment	the comment string, terminated with Hex:00

### RECORD TYPE 9: THE MEASUREMENT COMMENT RECORD

This record stores a comment. For the ICS standard it should be used to identify the patient and the type of the measurement.

Record type:	9	:integer
Blocksize:	N	:integer
Checksum:	$-(9 + \text{blocksize})$	:integer
Data block:	measurement comment	:string[N]

The following information should be stored:

In the field:

measurement comment a free-format string, terminated with Hex:00, that uniquely identifies patient and measurement for the originating center. This string would include not only a patient ID, e.g., from the local hospital information system, but also the relation of this measurement to a particular research study, e.g., "post-TURP investigation". For use in a specific multicenter study, this string can be specified more strictly.

### RECORD TYPE 14: THE ICS SIGNAL PROPERTY RECORD

This record type describes the properties of one stored signal, including its name, the unit of measurement, the zero and full-scale values, and the sample rate (which, if possible, should be identical for all signals).

Record Type:	14	:integer
Blocksize:	$1 + N + M + 3*2 + 4$	:integer
Checksum:	$-(14 + \text{blocksize})$	:integer
Data block:	signal id	:byte
	signal name	:string [N]
	unit	:string [M]
	binzero	:integer
	binsize	:word
	fullscale	:integer
	sample rate	:dword

The following information should be stored in these fields:

In the field:

signal id as specified in Section A9: Addendum  
signal name as specified in Section A9: Addendum

unit	as specified in Section A9: Addendum
binzero	the signal value represented in the data samples field of record type 18 by the word 0
binsize	the full scale binary sample value
fullscale	the signal value represented in the data samples field of record type 18 by the full scale binary value
sample rate	sample rate specification in samples/s (Hz). The sample rate is in 16.16 format: i.e., it is a long integer, where bits 16 through 31 specify the integer part of a fixed point number and bits 0 through 15 specify the decimal part.

The following values are recommended:

- sample rate at least 10 samples/s
- pressure signals bipolar up to +200 cm H<sub>2</sub>O with a minimum resolution of 0.5 cm H<sub>2</sub>O/bit
- flow rate signal unipolar up to 50 ml/s with a minimum resolution of 0.05 ml/s/bit

*Example* showing how sample rate is to be specified: Suppose the rate of sampling is 512.5 Hz. Then the value to be stored in the field sample rate is the 16.16 value Hex:02008000 = 512.5\*Hex: 10000.

*Example* showing how binzero, binsize, and fullscale are to be used: Suppose the abdominal pressure is sampled bipolarly with a 10-bit A/D converter. Then the binary value will range from 0 to 1023. Suppose the A/D converter is configured so that this range represents the pressure range from -50 to 200 cm H<sub>2</sub>O. Then the three values to be stored are:

```
binzero    =  -50
binsize    = 1023
fullscale  =  200
```

### RECORD TYPE 15: THE ICS PATIENT DATA RECORD

This record type contains the basic demographic data of the patient. It contains the actual text and numerical fields, plus an index table to facilitate location of a particular field. The index gives the relative position within the data block.

Record type:	15	:integer
Blocksize:	17*2 + N + M + O + P + Q + R + S + T + U + V + W + X + Y + 3 + 2	:integer
Checksum:	-(15 + blocksize)	:integer
Data block:	surname index	:word
	first name index	:word
	maiden name index	:word
	ID index	:word
	street index	:word
	houenumber index	:word
	city index	:word
	country index	:word
	postcode index	:word
	phone index	:word
	height index	:word
	weight index	:word
	sex index	:word
	comments index	:word

## 14 Griffiths et al.

birth day index	:word
birth month index	:word
birth year index	:word
surname	:string[N]
first name	:string[M]
maiden name	:string[O]
ID	:string[P]
street	:string[Q]
houzenumber	:string[R]
city	:string[S]
country	:string[T]
postcode	:string[U]
phone	:string[V]
height	:string[W]
weight	:string[X]
comments	:string[Y]
sex	:byte*
birth day	:byte
birth month	:byte
birth year	:word <sup>†</sup>

\*0 = male, 1 = female

<sup>†</sup>birth year is expressed in full, e.g., 1895, 1995, 2005

The first name field may contain a middle or other initial in addition to the first name.

As a default, an empty string (containing only Hex:00) may be stored in any of the string fields for which information is not available or may not be disclosed because of ethical considerations. In this case the length (N, M, etc.) of the string is 1.

The index fields contain the position (in bytes) of the corresponding text or numerical field, relative to the start of the data block. Thus the value to be stored in surname index is  $17*2$ ; the value to be stored in first name index is  $17*2 + N$ ; and so on. All the index fields must be completed.

### RECORD TYPE 16: THE ICS SOURCE RECORD

This record type specifies the origin of the measurement. It consists of an index table together with the actual text fields.

Record type:	16	:integer
Blocksize:	$10*2 + N + M + O + P + Q + R + S + T + U + V$	:integer
Checksum:	$-(16 + \text{blocksize})$	:integer
Data block:	clinic name index	:word
	investigator name index	:word
	street index	:word
	streetnumber index	:word
	city index	:word
	country index	:word
	postcode index	:word
	phone index	:word
	fax index	:word
	comments index	:word

clinic name	:string[N]
investigator name	:string[M]
street	:string[O]
streetnumber	:string[P]
city	:string[Q]
country	:string[R]
postcode	:string[S]
phone	:string[T]
fax	:string[U]
comments	:string[V]

The phone and fax fields should include area codes as well as local numbers.

An empty string may be stored in any of the string fields for which information is not available. In this case the length (N, M, etc.) of the string is 1.

The index fields contain the position (in bytes) of the corresponding text field, relative to the start of the data block. All the index fields must be completed.

**RECORD TYPE 17: ICS VOLUME RECORD**

This record contains the filling volume and the residual volume.

Record type:	17	:integer
Blocksize:	4	:integer
Checksum:	-(17 + blocksize)	:integer
Data block:	filling volume	:integer
	residual volume	:integer

The following information should be stored:

In the field:	
filling volume	the calculated volume in the bladder at the beginning of the pressure-flow study (in ml). -1 = default value
residual volume	the volume in the bladder at the end of the study, either calculated or measured directly (in ml); if the calculated volume is negative the value zero should be stored. -1 = default value

**RECORD TYPE 18: THE ICS SIGNAL VALUE RECORD**

The ICS signal value record contains the actual data samples for one of the urodynamic signals. It includes a record number, allowing a signal that is too long to fit into one record to be divided so as to span several records. If the signal spans more than one record, the number of the first record should be 1 and the following records should be numbered 2, 3, and so on. If there is only one record, its record number should be 0.

Record type:	18	:integer
Blocksize:	2 + N*2	:integer
Checksum:	-(18 + blocksize)	:integer
Data block:	signal id	:byte
	record number	:byte
	data samples	:word[N]

## 16 Griffiths et al.

The following information should be stored in these fields:

In the field:  
signal id as specified in Section A9: Addendum  
record number as described above  
data samples the binary samples themselves,  
1 word = 2 bytes = 1 sample,  
stored in the order low byte, high byte

### RECORD TYPE 19: THE ICS MEASUREMENT DESCRIPTION RECORD

This record indicates that the file is an ICS standard file and describes the measurement's start date and time, the number of signals stored, and the number of records in the file.

Record type:		:integer
Blocksize:	1 + N + 5 + 1 + 2	:integer
Checksum:	-(19 + blocksize)	:integer
Data block:	measurement type	:byte
	name	:string[N]
	start	:minute
		:hour
		:day
		:month
		:year
	number of signals	:byte
	number of records	:word

\*year is expressed modulo 100; e.g. 1995 = 95, 2005 = 05.

The following information should be stored in these fields:

In the field:  
measurement type the version number of this document, multiplied by 10 (e.g., for this version 71). Note that the version number may not contain hundredths or smaller decimal fractions  
name the string "ICS standard pressure-flow study"  
start the start time and date of the measurement (if any of these are unavailable, the default value Hex:00 should be stored in the corresponding position: minute, hour, day, month and/or year)  
number of signals the number of signals stored in this file (e.g., 3 if just intravesical pressure, abdominal pressure and flow rate are stored).  
number of records the total number of records stored in this file (including this measurement description record)

### A6. SIGNALS AND INFORMATION TO BE STORED: MINIMAL SPECIFICATION AND OPTIONAL EXTENSIONS

Minimally, the following 3 signals should be stored in 3 records of type 18:  
intravesical pressure  
abdominal pressure  
flow rate.

For each signal, certain associated information should be stored in a record of type 14.

Optionally, the following signals can also be stored in records of type 18, with associated information in records of type 14:

EMG envelope  
voided volume

The voided volume signal may be useful if this is the signal measured by the urodynamic system, and the flow signal is derived from it.

In addition to the signals, further information about the patient and the measurements should be stored in records of types 16, 17, and 19. Full demographic data for the patient can be stored in a record of type 15. Optionally, a free format comment can be stored in a record of type 9, and if detailed comments relating to events during the measurement are available, they can be stored in records of type 8.

## A7. TYPICAL FILE STRUCTURE

Of the various record types, some are mandatory and some are optional. The structure of the file describing a particular measurement varies according to what records are stored.

The order of the records within the file is arbitrary, in principle. However, it may be convenient to place records containing easily recognizable strings near the beginning, so that the file can quickly be identified if it is accidentally renamed or misplaced. In particular, record type 19 should preferably be the first one.

Thus a typical file structure might be:

RECORD TYPE 19: ICS MEASUREMENT DESCRIPTION RECORD (identifies file as an ICS standard file)

RECORD TYPE 16: ICS SOURCE RECORD (identifies originating clinic and investigator)

RECORD TYPE 9: MEASUREMENT COMMENT RECORD (identifies patient and type of measurement)

OPTIONALLY, RECORD TYPE 15: ICS PATIENT DATA RECORD  
(contains full demographic data for patient)

RECORD TYPE 14: ICS SIGNAL PROPERTY RECORD (describes properties of intravesical pressure signal)

RECORD TYPE 14: ICS SIGNAL PROPERTY RECORD (describes properties of abdominal pressure signal)

RECORD TYPE 14: ICS SIGNAL PROPERTY RECORD (describes properties of flow rate signal)

OPTIONALLY, further records of type 14 for voided volume and EMG envelope

RECORD TYPE 18: ICS SIGNAL VALUE RECORD (contains actual intravesical pressure data)

RECORD TYPE 18: ICS SIGNAL VALUE RECORD (contains actual abdominal pressure data)

RECORD TYPE 18: ICS SIGNAL VALUE RECORD (contains actual flow rate data)

OPTIONALLY, further records of type 18 for voided volume and EMG envelope

OPTIONALLY, RECORD TYPE 8: MEASUREMENT MARKER RECORD (identifies position and type of marker in data, the signal with which it is associated, and a comment)

OPTIONALLY, further records of type 8 with additional markers

RECORD TYPE 17: ICS VOLUME RECORD (contains filling volume prior to pressure-flow study and residual urine volume after voiding)

## A8. ACKNOWLEDGMENTS

This proposal is based partly on the ICS Measure File Format (ICSMFF) version 1.00 of 14.10.93, which was written by Michael Gondy Jensen, formerly of the Wiest Company and now of Andromeda Medical Systems, and also on an earlier proposal for a simple digital pressure-flow standard circulated by Ron van Mastrigt. Dantec Medical, Laborie Medical Technologies, Life-Tech, and Andromeda Medical Systems have agreed in principle to support a standard similar to the ICSMFF. The Life-Tech Company has agreed to support the approved ICS standard.

## A9. ADDENDUM: SIGNAL ID'S

In records of types 8, 14 and 18, a signal ID of the type byte is specified. Corresponding signal names and units of the type string are required in record type 14. The following signal ID values, signal names and units have been defined.

	Signal ID	Signal name	Unit
For intravesical pressure	1	“pves”	“cmH2O”
For abdominal pressure	2	“pabd”	“cmH2O”
For flow rate	3	“Q”	“ml/s”
For EMG envelope	4	“EMGenv”	“uV”
For voided volume	5	“Vvoid”	“ml”

The following are reserved for possible future use:

	Signal ID	Signal name	Unit
For $p_{det}$ acquired independently	6	“pdet”	“cmH2O”
For infused volume	7	“Vinf”	“ml”
For direct EMG signal	8	“EMG”	“uV”
For urethral pressure	9	“pura”	“cmH2O”
For urethral closure pressure	10	“pclos”	“cmH2O”

Additional signal ID's, signal names, and units may be introduced; proposed additions should be communicated to the ICS for registration.